HDMI Cameras – H Series

# INSTRUCTION MANUAL

<b>Model</b>
C-HP

v 1.0 2018

---

# Contents

C-HP

C-mount HDMI+USB CMOS Camera

Instruction Manual

- 1 C-HP Basic Characteristics
  - 1.1 USB Video Mode
  - 1.2 HDMI MODE.
- 2 Brief Introduction of C-HP UI and Its Functions
  - 2.1 The Camera Control Panel on the Left Side of the Video Window
  - 2.2 Icons and Functions of the Camera Control Toolbar at the Bottom of the Video Window
  - 2.3 The Measurement Toolbar on the Upper Side of the Video Window

---

# 1 C.HP Basic Characteristics

C-HP is a multiple interfaces (HDMI + USB2.0 + SD card) CMOS camera. HDMI + USB2.0 are used as the data transfer interface to HDMI display or computer.

For HDMI output, the **Camera Control Panel**+ **Measurement Toolbar** and **Camera Control Toolbar** are overlaid on the HDMI screen when the mouse move to the related region, in this case, the USB mouse can be used to set the camera, browse and compare the captured image, play the video and perform the measurement ital.

For **USB Video** camera mode, plug in the micro USB host cable to the camera **USB video** port and computer USB port, then the video stream can be transfer to computer with the advanced software ProView. With ProView, you can control the camera, process the video and image as OPTIKA's other USB series camera.

The C-HP's basic characteristics are as follows:

For HDMI output:

- All in 1(HDMI + USB + SD card) C-mount camera with Sony high sensitivity CMOS sensor;
- Simultaneous HDMI & USB output;
- Built-in mouse control;
- Built-in image capture & video record to SD card;
- Built-in camera control panel, including exposure(manual/auto)/gain, white balance(lockable), color adjustment, sharpness control;
- Built-in video and image measurement;
- Built-in toolbar including zoom, mirror, comparison, freeze, cross, browser functions;
- Built-in image & video browsing, display & play;
- Real time clock(RTC)

For USB Video output:

- Ultra-Fine color engine with perfect color reproduction capability(USB);
- With advanced video & image processing application ProView, which including professional image processing such as 2D measurement, HDR, image stitching, EDF(Extended Depth of Focus), image segmentation & count, image stacking, color composite and denoising (USB);
- Support standard UVC protocol for Windows/Linux/Mac(USB);

### C-HP - Microscope + Display

C-HP - SW UI for Mouse Control

Figure 1 - The Layout of C-HP Rear Cover

Before starting the camera, please connect the standard **C-HP** C-mount camera to the camera adapter and insert it into the microscope's 3<sup>rd</sup> photo tube which will relay microscope object's middle image to the camera sensor.

## 1.1 USB Video Mode

1. Plug **12V/1A** power cable into DC 12V slot to supply power for the camera. The **LED Indicator** will turn into red;
2. Press **ON/OFF Button** to start the camera and the **LED Indicator** will turn into blue;
3. Plug the USB cable which comes with the camera into **USB Video** to connect the video to the computer;
4. Open **ProView** software, start **C-HP** by clicking the camera model name listed in **Camera List**. For more details please refer to the **ProView** help manual;

## 1.2 HDMI MODE

1. Plug the HDMI cable into the **HDMI Port** to connect the **C-HP** camera to **HDMI** display;
2. Plug a USB mouse into **USB Mouse** to get control of the camera by using built-in software ProView;
3. Plug 12V/ 1A power adapter into **Power Interface** to supply power for the camera. The **LED Indicator** will turn into red;
4. Insert SD card into **SD Card Slot** for saving captured images and recorded videos;
5. Press **ON/OFF Button** to start the camera. The **LED Indicator** will turn into blue;
6. Move mouse cursor to the left side of the video window, a **Camera Control Panel** will appear. It includes **Manual/ Automatic Exposure**, **White Balance**, **Sharpness** and other functions, please refer to 2.1 for details;
7. Move mouse cursor to the bottom of the video window and a **Camera Control Toolbar** will appear. Operations like **Zoom In**, **Zoom Out**, **Flip**, **Freeze**, **Cross Line**, **Comparison** and etc. can be realized. Please refer to 2.2 for details;
8. Move mouse cursor to the upper side of the video window, a **Measurement Toolbar** with calibration and other measurement tools will appear, please refer to 2.3 for details; The measurement data can be output with **\*.CSV** format.

## 2 Brief Introduction of C-HP UI and Its Functions

The C-HP UI shown in Fig. 2 includes a **Camera Control Panel** on the left side of the video window, a **Measurement Toolbar** on the upper side of the video window, a **Camera Control Toolbar** on the bottom of the video window.



Figure 2 - The C-HP Camera Control UI

Notes	
1	When users move mouse cursor to the left side of the video window, the <b>Camera Control Panel</b> will pop up automatically;
2	When users move mouse cursor to the bottom of the video window, the <b>Camera Control Toolbar</b> will pop up automatically;
3	When user moves mouse cursor to the bottom of the video window, the <b>Camera Control Toolbar</b> will pop up automatically.
4	Move the mouse cursor to the upper side of the video window, a <b>Measurement Toolbar</b> will pop up for the calibration and measurement operations. When user left-clicks the <b>Float/Fixed</b> button  on the <b>Measurement Toolbar</b> , the <b>Measurement Toolbar</b> will be fixed. In this case the <b>Camera Control Panel</b> will not pop up automatically even if users move mouse cursor to left side of the video windows. Only when user left-clicks the  button on the <b>Measurement Toolbar</b> to exit from measuring procedure will they be able to do other operations on the <b>Camera Control Panel</b> , or <b>Synthesis Camera Control Toolbar</b> . During the measuring process, when a specific measuring object is selected an <b>Object Location &amp; Attributes Control Bar</b>  will appear for changing location and properties of the selected objects.

## 2.1 The Camera Control Panel on the Left Side of the Video Window

The **Camera Control Panel** controls the camera to achieve the best image quality according to the specific applications; It will pop up automatically when mouse cursor is moved to the left side of the video window (in measurement status, the **Camera Control Panel** will not pop up. Only when measurement process is terminated will the **Camera Control Panel** pop up by moving mouse cursor to the left side of the video window). Left-clicking  button to achieve **Display/ Auto Hide** switch of the **Camera Control Panel**;

Camera Control Panel	Function	Function Description
	<b>Snap</b>	<b>Capture</b> image from the current video window
	<b>Record</b>	<b>Record</b> video from the current video window
	<b>Auto Exposure</b>	When <b>Auto Exposure</b> is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	<b>Target</b>	Available when <b>Auto Exposure</b> is checked. Slide to left or right to adjust <b>Target</b> according to the current video brightness to achieve proper brightness value
	<b>Exposure Time</b>	Available when <b>Auto Exposure</b> is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	<b>Gain</b>	Adjust <b>Gain</b> to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	<b>Red</b>	Slide to left or right to decrease or increase the proportion of <b>Red</b> in RGB on video
	<b>Green</b>	<b>Green</b> is base for reference and cannot be adjusted
	<b>Blue</b>	Slide to left or right to decrease or increase the proportion of <b>Blue</b> in RGB on the video
	<b>White Balance</b>	<b>Auto</b> : <b>White Balance</b> adjustment according to the window video; <b>Manual</b> : Slide the <b>Red</b> or <b>Blue</b> to manually set the video <b>White Balance</b> ; <b>ROI</b> : Set the <b>White Balance</b> according to the <b>ROI</b> . The <b>ROI</b> can be resized and moved;
	<b>Sharpness</b>	Adjust <b>Sharpness</b> level of the video window
	<b>Saturation</b>	Adjust <b>Saturation</b> level of the video window
	<b>Gamma</b>	Adjust <b>Gamma</b> level of the video. Slide to the right side to increase gamma and to the left to decrease gamma.
	<b>Contrast</b>	Adjust <b>Contrast</b> level of the video. Slide to the right side to increase contrast and to the left to decrease contrast.
	<b>DC</b>	For <b>DC</b> illumination, there will be no fluctuation in light source so no need for compensating light flickering
	<b>AC(50HZ)</b>	Check <b>AC(50HZ)</b> to eliminate flickering “strap” caused by 50Hz illumination
<b>AC(60HZ)</b>	Check <b>AC(60HZ)</b> to eliminate flickering “strap” caused by 60Hz illumination	
<b>Default</b>	Set all the settings in the <b>Camera Control Panel</b> to default values	

## 2.2 Icons and Functions of the Camera Control Toolbar at the Bottom of the Video Window

Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Reset Zoom to Original 1X		Horizontal Flip
	Vertical Flip		Video Freeze
	Compare Image with the current video		Display Cross Line
	Browse Images and Videos in the SD Card		Settings
	Check the Version of ProView		

The Setting function is relatively more complicated than the other functions. Here are more info about it:


Figure 3 - Comprehensive Measurement Settings Page

**Global** : Used for setting digits behind the decimal point for measurement results;

**Calibration LineWidth**: Used for defining width of the lines for calibration;

**Color**: Used for defining color of the lines for calibration;

**EndPoint** Type: Used for defining shape of the endpoints of lines for calibration: Null means no endpoints, rectangle means rectangle type of endpoints. It makes alignment more easily;

**Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve**:

Left-click the besides the measuring patterns mentioned above will unfold the corresponding attribute settings to set the individual property of the measuring objects.


Figure 4 - Comprehensive Magnification Calibration Management Settings Page

**Name:** Names such as 10X, 40X, 100X are based on magnification of the microscopes. For continuous zoom microscopes, ensure that the selected magnification coincides with the scale alignment line on the microscope zoom knob;

**Resolution:** Pixels per meter. Devices like microscopes have high resolution value;

**Clear All:** Click the **Clear All** button will clear the calibrated magnifications and resolutions;

**Delete:** Click **Delete** to delete the selected item for specific resolution;

Figure 5 - Comprehensive Setting of SD Card Setting Page

**Current File System:** The maximum file **FAT32** can store is of 4G Bytes; for **EXFAT**, it's 2048G Bytes. Suggest converting **FAT32** file into **EXFAT** format on a PC; **Unknown Status:** SD card not detected or the file system is not identified;

Figure 6 - C-HP Comprehensive Setting of Language Selection Setting Page

- English :** Set language of the whole software into English;
- Simplified Chinese :** Set language of the whole software into Simplified Chinese;
- Traditional Chinese :** Set language of the whole software into Traditional Chinese;
- Korean:** Set language of the whole software into Korean;
- Thailand:** Set language of the whole software into Thailand;

Figure 7 - The Real Time Clock Setting

- Year:** Current Year
- Month:** Current Month
- Day:** Current Day
- Hour:** Current Hour
- Minute:** Current Minute
- Second:** Current Second

## 2.3 The Measurement Toolbar on the Upper Side of the Video Window

The **Measurement Toolbar** will pop up when moving mouse cursor to any place near the upper side of the video window. Here are the introduction of the various functions on the **Measurement Toolbar**:

Figure 8 - The Measurement Toolbar Button on the Upper Side of the Video window

Icon	Function
	Float/ Fix switch of the <b>Measurement Toolbar</b>
<input checked="" type="checkbox"/> Visible	Define measuring object in <b>Show up/ Hide</b> mode
Nanometer (nm)	Select the desired <b>Measurement Unit</b>
4X	Choose the same <b>Magnification</b> as the microscope to ensure accuracy of measurement result when measurement unit is not in Pixel unite
	<b>Object Select</b>
	<b>Point</b>
	<b>Angle</b>
	<b>Arbitrary Line</b>
	<b>Parallel</b>
	<b>Horizontal Line</b>
	<b>Vertical Line</b>
	<b>Rectangle</b>
	<b>Circle</b>
	<b>Ellipse</b>
	<b>Annulus</b>
	<b>Two Circles and Center Distance</b>
	<b>Arc</b>
	<b>Polygon</b>
	<b>Curve</b>
	<b>Make Calibration</b> to determine the corresponding relation between magnification and resolution, this will establish the corresponding relationship between measurement unit and the sensor pixel size. Calibration needs to be done with the help of a micrometer. For detailed steps of carrying out calibration please refer to ProView help manual.
	Export the measurement information to <b>CSV file(*.csv)</b>
	<b>Delete All the Measurement Objects</b>
	<b>Setting</b>

	<p>Exit from Current Measurement Mode</p>
	<p>When the measurement ends, left-click on a single measuring object and the <b>Object Location &amp; Properties Control Bar</b> will show up. The icons on the control bar mean <b>Move Left</b>, <b>Move Right</b>, <b>Move Up</b>, <b>Move Down</b>, <b>Color Adjustment</b> and <b>Delete</b>.</p>

Note:

1) When user left-clicks **Display/Hide** button  on the **Measurement Toolbar**, the **Measurement Toolbar** will be fixed. In this case the **Camera Control Panel** will not pop up automatically even if moving mouse cursor to the left side of the video window. Only when users left click the  button on the **Measurement Toolbar** to exit from the measurement mode will they be able to doing other operations in the **Camera Control Panel**, the **Auto Focus Control Panel** or the **Synthesis Camera Control Toolbar**.

2) When a specific measuring object is selected during the measuring process, the **Object Location & Attributes Control Bar**

 will appear for changing the object location and properties of the selected objects

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**OPTIKA® S.r.l.**